## INTENDED USE

Feline Herpes virus (FHV) is a common causal agent of upper respiratory disease and eye inflammation in felines. Asan Easy Test<sup>®</sup> FHV is a chromatographic immunoassay kit for rapid and qualitative detection of herpes virus from in feline eye mucus and nasal discharge. Antigens of herpes virus in the specimens are allowed to react with the anti-FHV monoclonal antibody-coupled gold conjugate followed by reaction with another anti-FHV monoclonal antibodies immobilized in the test lines. The result is visible within 10 minutes in the form of a red or purple line on the test line (T) area of the device. The solution continues to migrate to encounter a control reagent that binds a control conjugate, thereby producing a second red line on the control line (C) area.

## MATERIALS PROVIDED

Asan Easy Test<sup>®</sup> FHV contains the following components (10Test/Kit)

1. Test device individually foil-pouched with a desiccant : 10units
2. Assay solution in tube (0.5ml/tube) : 10units
3. Sterilized swab : 10units
4. Disposable dropper : 10units
5. Instruction manual for use : 1sheet

## PRECAUTIONS


1. The presence of humidity may decrease the stability of the reagents.
2. Do not use the kit after the expiration date and do not freeze.
3. For *in vitro* diagnostic use only.
4. Wear protective gloves while handling samples and wash hands thoroughly after the test.
5. Dispose all the specimens and kits properly after test, in accordance with GLP.
6. Never use reagents from another kit.
7. Discard the assay solution if it contaminated with bacteria or mold.

## SPECIMEN COLLECTION AND STORAGE

1. The test should be performed using feline eye mucus infected with conjunctival fornix or nasal discharge.
2. Specimen to be tested should be obtained and handled by standard methods for the collection of body secretions.
3. After collecting the specimen using swab, all the specimens should be immediately extracted and tested. If necessary, they may be stored at 2 to 8°C for upto 2 days or at -20°C for longer periods


## TEST PROCEDURE

1. Place all specimens, test devices, and allow them to room temperature prior to testing (15~30 min).
2. Take a portion of specimen using the sterilized swab.
3. Insert the swab sample into the test tube, then swirl it at least 6 times while pressing the head against the bottom and side of the test tube.
4. Roll the swab head against the inside of the tube as you remove it. Dispose of the used swab in accordance with your biohazard waste disposal protocol.
5. Load 3~4 drops of the solution taken by the disposable dropper into the sample well(S).
6. Interpret the test result within 10 minutes. Do not read the result after 10 minutes.

## INTERPRETATION OF THE TEST

1. Negative result : ONLY one band in the control line (C).
2. Positive result : Two bands are appeared in the test line (T) and control line (C).
3. Invalid result : If at 10 minutes, the red color band does not appear in the control line (C), even if any shade of a pink-to-red test line (T) appears, the result is considered invalid. If the test is invalid, a new test should be performed with a new sample and a new test device.

## STORAGE & EXPIRATION

1. Asan Easy Test<sup>®</sup> FHV kit should be stored between 1 to 30°C (33.8~86°F).
2. Expiration date of this kit is 24 months after its manufacture date.

## PERFORMANCE CHARACTERISTIC

Asan Easy Test<sup>®</sup> FHV is designed for primary screening test of feline herpes virus antigen. This kit can provide fast and easy way to get a result, but do not completely exclude the possibility of false positive or false negative result caused by various factors. So, refer to the result of this kit, please make a final decision with clinical manifestation, other test results, and veterinarian's view, collectively.

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