For Veterinary Use Only

Asan Easy Test[®] TgAb

Toxoplasma gondii antibody (Tg Ab) test in serum, plasma and whole blood | 林特睿检™犬、猫弓形虫抗体快速诊断试纸
Immunochromatography

EXPLANATION OF THE TEST

The obligate intracellular protozoan *Toxoplasma gondii* is ubiquitously distributed and infects a wide range of warm-blooded host, including up to 30% of human population worldwide. Infection of immunocompetent humans is most often asymptomatic, but leads to lifelong persistence of the parasite. *T. gondii* may lead to life-threatening disease in fetuses or new borns from primarily infected mothers or after reactivation of dormant parasites in immunocompromised patients such as AIDS or other immunosuppressive therapy.

Asan Easy Test[®] TgAb is a chromatographic immunoassay kit for rapid and qualitative detection of all antibodies (IgG, IgM, IgA) against *T. gondii* in serum, plasma, and whole blood from feline or canine. The nitrocellulose membrane of the kit is immobilized with *T. gondii*-specific antigens in the test line and with anti-nus-tag in the control line. And also, *T. gondii*-specific antigens and nus-tag are conjugated to the colloidal gold particles. This conjugate is placed on a polyester or glass fiber pad as conjugate pad. When the sample is dropped into the sample well on the device, the solubilized conjugate migrates with the sample by passive diffusion and both the conjugate and sample come into contact with the antigens that immobilized onto the nitrocellulose. If the sample contains antibodies against *T. gondii*, the result is visible as red line within 10 minutes in the test line on the membrane. The solution continues to migrate to encounter a control reagent that binds a control conjugate, thereby producing another red control line.

MATERIALS PROVIDED

- Asan Easy Test[®] TgAb kit contains the following components (10 Test/Kit):
1. Test device individually foil-pouched with a desiccant: 10 units
 2. Assay solution (2 mL/bottle): 1 unit
 3. Capillary pipette: 10 units
 4. EDTA tube: 10 units
 5. Instruction manual for use: 1 sheet

PRECAUTIONS

1. The presence of humidity may decrease the stability of the reagents.
2. Do not use the kit after the expiration date and do not freeze.
3. For veterinary use only.
4. Wear protective gloves while handling samples and wash hands thoroughly after the test.
5. Dispose all the specimens and kits properly after test, in accordance with GLP.
6. Never use reagents from another kit.
7. Discard the Assay solution if it is contaminated with bacteria or mold.

SPECIMEN COLLECTION AND STORAGE


1. Specimen to be tested should be obtained and handled by standard methods for their collections.

- Serum:** Allow the blood to clot, then centrifuge to separate the serum.
- Plasma:** Collect the whole blood into the tube contained anticoagulants such as heparin, citrate, or EDTA. Centrifuge the blood and separate the plasma.
- Whole blood:** Collect an anticoagulated blood sample in EDTA, heparin or citrate using standard clinical laboratory procedures. Anticoagulated whole blood samples should be tested within 24 hours of drawing. If delays are expected between sampling and testing, the sample should be stored either on ice or refrigerated (2-8°C), but should not be frozen.

2. All specimens should be tested as soon as they are prepared. If necessary, they may be stored at 2-8°C for up to 24 hours or at -20°C for longer periods.

TEST PROCEDURE

1. Place all specimens, test devices and assay solution and allow them to room temperature prior to testing (15-30 min).
2. Prepare the test device as you need, and then mark the dog's ID onto the device. Please perform the test immediately after removing the device from foil pouch.
3. With 10 μ L capillary pipette, add 10 μ L of specimen drawn to black line into the sample well (S).
4. Add 3-4 (approx. 100 μ L) drops of Assay solution into the sample well (S).
5. Interpret the testing result within 10 minutes. Do not read the result after 10 minutes.

INTERPRETATION OF THE RESULTS


1. Negative result: ONLY one band in the control line (C).

2. Positive result

Two bands are appeared in the test line (T) and control line (C).

3. Invalid result:

If at 10 minutes, the red color band does not appear in the control line (C), even if any shade of a pink-to-red test line (T) appears, the result is considered invalid. If the test is invalid, a new test should be performed with a new sample and a new test device.

STORAGE & EXPIRATION

1. Asan Easy Test[®] TgAb kit should be stored between 1 to 30°C (33.8-86°F).
2. Expiration date of this kit is 24 months after its manufacture date.

LIMITATIONS OF THE TEST

Asan Easy Test[®] TgAb is designed for primary screening test of antibodies against *T. gondii*. This kit can provide fast and easy way to get a result, but do not completely exclude the possibility of false positive or false negative result caused by various factors. So, refer to the result of this kit, please make a final decision with clinical manifestation, other test results, and veterinarian's view, collectively.

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