For Veterinary Use Only

Asan Easy Test[®] Heartworm

Canine Heartworm Virus Antigen (CHW Ag) Test | 林特睿检[™]犬心丝虫抗原快速诊断试纸

Immunochromatography

EXPLANATION OF THE TEST

Asan Easy Test[®] Heartworm is a chromatographic immunoassay kit for rapid and qualitative detection of dirofilariasis antigen in canine whole blood, serum or plasma. A nitrocellulose membrane is immobilized with antibody against the D. immitis antigen. Another anti-D. immitis antibody is conjugated to colloidal gold particles. This conjugate is placed on a polyester pad as a conjugate pad. This test is aimed to the detection of D. immitis antigen in the blood. When the device is loaded with assay solution containing blood specimens, the solubilized conjugate migrates with the specimens by passive diffusion and both the conjugate and sample come into contact with the anti-D. immitis antibody that it adsorbed onto the nitrocellulose. The result is visible within 10 minutes in the form of a red or purple line on the test line (T) region of the device. The solution continues to migrate to encounter a control reagent that binds a control conjugate, thereby producing a second red line on the control line (C) area.

MATERIALS PROVIDED

- Asan Easy Test[®] Heartworm contains the following components : (20Tests/Kit)
 1. Test devices individually foil-pouched with a desiccant : 20units
 2. Assay solution (5ml/bottle) : 1unit
 3. EDTA Tube : 20units
 4. Disposable dropper : 20units
 5. Instruction manual for use : 1sheet

PRECAUTIONS


1. The presence of humidity may decrease the stability of the reagents. Thus, please carry out the test immediately after removing the device from the foil pouch.
2. Do not use the kit after the expiration date and do not freeze the kit.
3. For veterinary use only.
4. Wear protective gloves while handling samples and wash hands thoroughly after the test.
5. Dispose all the specimens and kits properly after test, in accordance with GLP.
6. Never use reagents from another kit.
7. Discard the Assay solution if it is contaminated with bacteria or mold.

SPECIMEN COLLECTION AND STORAGE

1. Specimen to be tested should be obtained and handled by standard methods for the collection of blood from canine.
 - A. Serum : Allow the blood to clot, then centrifuge to separate the serum.
 - B. Plasma : Collect the whole blood into the tube contained anticoagulants such as heparin, citrate, or EDTA. Centrifuge the blood and separate the plasma.
 - C. Whole blood : Collect an anticoagulated blood sample in EDTA, heparin or citrate using standard clinical laboratory procedures. Anticoagulated whole blood samples should be tested within 24 hours of drawing. If delays are expected between sampling and testing, the sample should be stored refrigerated (2~8°C), but should not be frozen.
2. Serum and plasma should be tested as soon as they are prepared. If necessary, they may be stored at 2~8°C for up to 24 hours or at -20°C for longer periods.


TEST PROCEDURE

1. Place all specimens, test devices, and allow them to room temperature prior to testing (15~30 min).
2. Remove the test device from the sealed pouch, and place it on a clean and flat place.
3. Take the blood sample (whole blood, serum or plasma) and load 1 drop (approx. 30µl) of the sample into sample well (S) of the device.
4. Hold the assay solution bottle vertically, add 3 drops of the solution into sample well (S) of the device.
5. Within 10 minutes, interpret the test result(s). Do not read the result after 10 minutes. Strong positive results may appear sooner.

INTERPRETATION OF THE RESULTS

1. Negative result : ONLY one band in the control line (C).
2. Positive result : Two bands are appeared in the test line (T) and control line(C).
3. Invalid result : If at 10 minutes, the red color band does not appear in the control line, even if any shade of a pink-to-red test line appears, the result is considered invalid. If the test is invalid, a new test should be performed with a new sample and a new test device.

STORAGE & EXPIRATION

1. Asan Easy Test[®] Heartworm kit should be stored between 2 to 30°C (36-86°F).
2. Expiration date of this kit is 24 months after its manufacture date.

LIMITATIONS OF THE TEST

Asan Easy Test[®] Heartworm is designed for primary screening test of canine heartworm antigen (Dirofilaria immitis antigen). This kit can provide fast and easy way to get a result, but do not completely exclude the possibility of false positive or false negative result caused by various factors. So, refer to the result of this kit, please make a final decision with clinical manifestation, other test results, and veterinarian's view, collectively.

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ASAN PHARMACEUTICAL CO., LTD