For Veterinary Use Only

Asan EasyTest[®] FeLV

Feline Leukemia Virus Antigen (FeLV Ag) Test | 林特睿检™猫白血病病毒抗原快速诊断试纸
Immunochromatography

EXPLANATION OF THE TEST

FeLV, or Feline leukemia virus, is a contagious, viral disease of cats. In addition to causing leukemia, it has been associated with various other types of cancer, anemia, and immune suppression leading to increased susceptibility to various infectious diseases. Although cats may clear initial infection, there is no cure for persistent infection and it is ultimately fatal. The virus is spread through direct contact with an FeLV positive cat's body fluids or saliva. It is disputed if the virus can live long enough outside of its host to spread itself from food or water bowls, or from a litter box. In any event, all cats owners should take precautions to protect their cats from the virus. They may live for as little as a few weeks, but probably not any longer than 3 years after their infection. But, even though they do not exhibit any signs, they will still transmit the disease to other cats - therefore the importance of a periodic test (or, even better, vaccination). Kittens of infected cats will most likely test positive for the virus, but signs of disease may not appear until later in their life. Asan EasyTest FeLV is a chromatographic immunoassay kit for rapid and qualitative detection of FeLV antigen in feline whole blood, plasma or serum. The nitrocellulose membrane of the kit is immobilized with FeLV-specific antibody in test line and with anti-mouse antibodies in the control line. And also, FeLV-specific antibody is conjugated to the colloidal gold particles. This conjugate is placed on a polyester or glass fiber as conjugate pad. When the sample is dropped into the sample well on the device, the solubilized conjugate migrates with the sample by passive diffusion and both the conjugate and sample come into contact with the antibody that immobilized onto the nitrocellulose. If the sample contains FeLV antigens, the result is visible as red line within ~10 minutes in the test line on the membrane. The solution continues to migrate to encounter a control reagent that binds a control conjugate, thereby producing another red control line.

MATERIALS PROVIDED

- Asan Easy Test FeLV kit contains the following components (10 Test/Kit):
1. Test device individually foil-pouched with a desiccant : 10units
 2. Assay solution in bottle (3ml/vial) : 1vial
 3. Capillary tube : 10units
 4. EDTA Tube : 10units
 5. Instruction manual for use : 1sheet

PRECAUTIONS


1. The presence of humidity may decrease the stability of the reagents.
2. Do not use the kit after the expiration date and do not freeze.
3. For veterinary use only.
4. Wear protective gloves while handling samples and wash hands thoroughly after the test.
5. Dispose all the specimens and kits properly after test, in accordance with GLP.
6. Never use reagents from another kit.
7. Discard the Assay solution if it is contaminated with bacteria or mold

SPECIMEN COLLECTION AND STORAGE

1. Specimen to be tested should be obtained and handled by standard methods for their collections.
 - A. Serum: Allow the blood to clot, then centrifuge to separate the serum.
 - B. Plasma: Collect the whole blood into the tube contained anticoagulants such as heparin, citrate, or EDTA. Centrifuge the blood and separate the plasma.
 - C. Whole blood: Collect an anticoagulated blood sample in EDTA, heparin or citrate using standard clinical laboratory procedures. Anticoagulated whole blood samples should be tested within 24 hours of drawing. If delays are expected between sampling and testing, the sample should be stored either on ice or refrigerated (2-8°C), but should not be frozen.
2. All specimens should be tested as soon as they are prepared. If necessary, they may be stored at 2-8°C.


TEST PROCEDURE

1. Place all specimens, test devices and Assay solution and allow them to room temperature prior to testing (15-30 min).
2. Prepare the test device as you need, and then mark the patients, ID onto the device. Please perform the test immediately after removing the device from foil pouch.
3. Load 1 drop (10 µl) of specimen into the sample well (S).
4. Add 4 drops of Assay Solution into the sample well (S).
5. Interpret the testing result within 10 minutes. Do not read the result after 10 minutes.

INTERPRETATION OF THE RESULTS

1. Negative result : ONLY one band in the control line (C).
2. Positive result : Two bands are appeared in the test line (T) and control line (C).
3. Invalid result : If at 10 minutes, the red color band does not appear in the control line (C), even if any shade of a pink-to-red test line (T) appears, the result is considered invalid. If the test is invalid, a new test should be performed with a new sample and a new test device.

STORAGE & EXPIRATION

1. Asan EasyTest FeLV kit should be stored between 1 to 30°C (33.8-86°F).
2. Expiration date of this kit is 24 months after its manufacture date.

LIMITATIONS OF THE TEST

Asan EasyTest FeLV is designed for primary screening test of feline leukemia virus infection. This kit can provide fast and easy way to get a result, but do not completely exclude the possibility of false positive or false negative result caused by various factors. So, refer to the result of this kit, please make a final decision with clinical manifestation, other test results, and veterinarian's view, collectively.

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ASAN PHARMACEUTICAL CO., LTD