For Veterinary Use Only

Asan Easy Test[®] Distemper

Canine Distemper Virus Antigen (CDV Ag) Test

Immunochromatography

EXPLANATION OF THE TEST

Asan Easy Test[®] Distemper is a chromatographic immunoassay kit for rapid and qualitative detection of body secretions of dog. A nitrocellulose membrane is immobilized with monoclonal antibody against the CDV antigen. Another anti-CDV Ag-monoclonal antibody is conjugated to colloidal gold particles. This conjugate is placed on a polyester pad as a conjugate pad. This test is aimed to the detection of CDV antigen either in canine body secretions or in virus-cultured supernatant. When the device is loaded with assay buffer containing blood specimen, the solubilized conjugate migrates with the sample by passive diffusion and both the conjugate and specimen come into contact with the anti-CDV Ag antibody that it adsorbed onto the nitrocellulose. The result is visible within 10 minutes in the form of a red or purple line on the test line (T) area of the device. The solution continues to migrate to encounter a control reagent that binds a control conjugate, thereby producing a second red line on the control line (C) area.

MATERIALS PROVIDED

Asan Easy Test[®] Distemper contains the following components

1. Test device individually foil-pouched with a desiccant
2. Assay solution in tube (0.5ml/tube)
3. Disposable dropper
4. Sterilized swab
5. Instruction manual for use

PRECAUTIONS

1. The presence of humidity may decrease the stability of the reagents. Thus, please carry out the test immediately after removing the device from the foil pouch.
2. Do not use the kit after the expiration date and do not freeze the kit.
3. For veterinary use only.
4. Wear protective gloves while handling samples and wash hands thoroughly after the test.
5. Dispose all the specimens and kits properly after test, in accordance with GLP.
6. Never use reagents from another kit.
7. Discard the Assay solution if it is contaminated with bacteria or mold.


SPECIMEN COLLECTION AND STORAGE

1. The test should be performed using canine secretion of eye mucus (in the conjunctiva) or nasal discharge.
2. Specimen to be tested should be obtained and handled by standard methods for the collection of body secretions.
3. After collecting the specimen using swab, all the specimens should be immediately extracted and tested. If necessary, they may be stored at 2~8°C for up to 48 hours or at -20°C for longer periods.

TEST PROCEDURE


1. Place all specimens, test devices, and allow them to room temperature prior to testing (15~30 min).
2. Remove the test device from the sealed pouch, and place it on a clean and plain place.
3. Take a portion of secretion of eye mucus or nasal discharge, using a sterilized swab.
4. Insert and mix the swab with assay solution in a tube until the specimens has been dissolved. Then, discard the swab and repeat procedure for each sample.
5. Add 3~4 drops (approx. 100 μ l) of the solution taken by the disposable dropper into the sample well.

6. Interpret the test result(s) within 5~10 minutes. Do not read the result after 10 minutes.

INTERPRETATION OF THE RESULTS

1. Negative result : ONLY one band in the control line (C).
2. Positive result : Two bands are appeared in the test line (T) and control line (C).
3. Invalid result : If at 10 minutes, the red color band does not appear in the control line, even if any shade of a pink-to-red test line appears, the result is considered invalid. If the test is invalid, a new test should be performed with a new sample and a new test device.

STORAGE & EXPIRATION

1. Asan Easy Test[®] Distemper kit should be stored between 2 to 30°C (35.6-86 °F).
2. Expiration date of this kit is 24 months after its manufacture date.

LIMITATION OF THE TEST

Asan Easy Test[®] Distemper is designed for primary screening test of canine distemper virus antigen. This kit can provide fast and easy way to get a result, but do not completely exclude the possibility of false positive or false negative result caused by various factors. So, refer to the result of this kit, please make a final decision with clinical manifestation, other test results, and veterinarian's view, collectively.

Manufactured & Sold by
ASAN PHARMACEUTICAL CO., LTD
 Factory1: 163, Yeongcheon-ro, Hwaseong-si, Gyeonggi-do 18462, Korea
 Factory2: 122-26, Gieopdanji-ro, Gongdo-eup, Anseong-si, Gyeonggi-do, 17551, Korea

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