## EXPLANATION OF THE TEST

Asan Easy Test® CCV/CPV is a chromatographic immunoassay kit for rapid and qualitative detection of canine coronavirus (CCV) and canine parvovirus (CPV) infection from feces specimen of dog. A nitrocellulose membrane is immobilized with monoclonal antibodies against the CCV and CPV antigens. Another anti-CCV Ag-monoclonal and anti-CPV Ag-monoclonal antibodies are conjugated to colloidal gold particles. These conjugates are placed on a polyester pad as a conjugate pad. This test is aimed to the detection of CCV and CPV antigens either in stool (feces) or in virus-cultured supernatant. When the device is loaded with assay buffer containing specimens, the solubilized conjugate migrates with the sample by passive diffusion and both the conjugate and sample come into contact with the anti-CCV and anti-CPV Ag antibodies that it adsorbed onto the nitrocellulose. The result is visible within 10 minutes in the form of a red or purple line on the test line (T1; CCV Ag, T2; CPV Ag) regions of the device. The solution continues to migrate to encounter a control reagent that binds a control conjugate, thereby producing a second red line on the control line (C) area.

## MATERIALS PROVIDED

Asan Easy Test® CCV/CPV contains the following components(10Test/Kit) :

1. Test devices individually foil-pouched with a desiccant : 10units
2. Assay solution in tube(1 mL/tube) : 10units
3. Disposable dropper : 10units
4. Sterilized swab : 10units
5. Instruction manual for use : 1sheet

## PRECAUTIONS


1. The presence of humidity may decrease the stability of the reagents.
2. Do not use the kit after the expiration date and do not freeze.
3. For in vitro diagnostic use only
4. Wear protective gloves while handling samples and wash hands thoroughly after the test
5. Dispose all the specimens and kits properly after test, in accordance with GLP.
6. Never use reagents from another kit.
7. Discard the Assay solution if it contaminated with bacteria or mold.

## SPECIMEN COLLECTION AND STORAGE

Specimen of this is canine fecal material(stool). And, it should be obtained and handled by methods for the sample collection. The feces specimens must be tested as soon early they are collected. If necessary, they may be stored at 2~8 °C for up to 48 hours. If longer storage is required, the sample may be frozen at -20 °C .


## TEST PROCEDURE

1. Place all specimens, test devices, and allow them to room temperature prior to testing (15~30 min).
2. Remove the test device from the sealed pouch, and place it on a clean and plain place.
3. Collect the feces specimens at 4 different points or directly take the specimens from the canine colon.
4. Insert and mix the swab with assay solution in a tube until the specimens has been dissolved. Then, discard the swab and repeat procedure for each sample.
5. Add 4 drops (approx. 120 μl) of the solution taken by the disposable dropper into the sample well.
6. Interpret the test result(s) within 10~15 minutes. Do not read after 15 minutes.

## INTERPRETATION OF THE RESULTS

1. Negative result : ONLY one band in the control line (C).
2. Positive result : Two bands are appeared in the test line (T) and control line(C).
3. Invalid result: If at 15 minutes, the red color band does not appear in the control line(C) , even if any shade of a pink-to-red test line(T) appears, the result is considered invalid. If the test is invalid, a new test should be performed with a new sample and a new test device.

## STORAGE & EXPIRATION

1. Asan Easy Test® CCV/CPV kit should be stored between 4 to 30 °C (39.2~86 °F) 4°C
2. Expiration date of this kit is 24 months after its manufacture date.

## LIMITATION OF THE TEST

Asan Easy Test® CCV/CPV is designed for primary screening test of canine coronavirus and canine parvovirus antigen. This kit can provide fast and easy way to get a result, but do not completely exclude the possibility of false positive or false negative result caused by various factors. So, refer to the result of this kit, please make a final decision with clinical manifestation, other test results, and veterinarian's view, collectively.

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